

Product datasheet for **RC208833A1V**

Human HOXC9 (NM_006897) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HOXC9 (NM_006897) AAV Particle
Tag:	Myc-DDK
Symbol:	HOXC9
Synonyms:	HOX3; HOX3B
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208833 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCGGCGACGGGGCCCATCAGTAACTATTACGTGGACTCGCTCATCTCTCACGACAATGAAGACCTCC
TAGCGTCCAGGTTTCCGGCCACCGGGGCTCATCCCGCCGCCAGACCCAGCGGTTTGGTGCCGGACTG
TAGCGATTTTCCGTCTGTAGCTTCGCGCCCAAGCCGGCAGTGTTTCAGCACGTCGTGGCGCCCCGTGCC
TCTCAGTCGTCCGTGGTATATCACCCGTACGGCCCCAGCCCCACCTCGGCGCCGACACGCGCTACATGC
GGACTTGGCTCGAGCCGCTGTCCGGCGCCGTCTCCTTCCCCAGCTTCCCGGCCGGGGCCGCTCACTACGC
CCTCAAGCCGGACGCCTACCCCGGGCGCGCGGACTGCGGCCAGGGGAGGGCCGAGCTACCCGGAC
TACATGTACGGCTCGCCCGGGGAGCTGCGCGACCGCGCCCGCAGACACTGCCCTCGCCCGAGGCGGACG
CGCTCGCCGGCAGCAAGCACAAGAGGAGAAGGCTGACCTGGACCCAGCAACCCCGTGGCCAACTGGAT
TCACGCCCCGCTCCACGAGGAAGAAGCGCTGCCCTACACCAAGTACCAGACGCTGGAAGTGGAGAAGGAG
TTTCTCTTCAATATGTATTTAACAGGGACCGTCGGTATGAGGTGGCCCGGGTTCTCAATCTCACCGAGC
GGCAGGTCAAAATCTGGTTTCAGAATCGAAGGATGAAGATGAAAAAGATGAATAAAGAGAAAACCGACAA
GGAGCAGTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208833 protein sequence Red=Cloning site Green=Tags(s) MSATGPISNYYVDSLISHDNEDLLASRFPATGAHPAAARPSGLVPDCSDFPSCSFAPKPAVFSTSWAPVP SQSSVVYHPYGPQPHLGADTRYMRTWLEPLSGAVSFPSFPAGGRHYALKPDAYPGRRADCGPGEGRSYPD YMYGSPGELRDRAPQTLPSPEADALAGSKHKEEKADLDPSNPVANWIHARSTRKKRCPYTKYQTLELEKE FLFNMYLTRDRRYEVARVLNLTERRQVKIWFQNRMRMKMKMNKEKTDKEQS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006897
ORF Size:	780 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006897.1
RefSeq Size:	1640 bp
RefSeq ORF:	783 bp
Locus ID:	3225
UniProt ID:	P31274
Cytogenetics:	12q13.13
MW:	29.2 kDa