

Product datasheet for **RC208846L3V**

NXPH4 (NM_007224) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NXPH4 (NM_007224) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NXPH4
Synonyms:	NPH4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007224
ORF Size:	924 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208846).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007224.1 , NP_009155.1
RefSeq Size:	1904 bp
RefSeq ORF:	927 bp
Locus ID:	11247
UniProt ID:	O95158
Cytogenetics:	12q13.3
Protein Families:	Secreted Protein
MW:	32.9 kDa



[View online »](#)

Gene Summary:

May be signaling molecules that resemble neuropeptides and that act by binding to alpha-neurexins and possibly other receptors.[UniProtKB/Swiss-Prot Function]