

Product datasheet for **RC208858A1V**

Human SPANX (SPANXC) (NM_022661) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SPANX (SPANXC) (NM_022661) AAV Particle
Tag: Myc-DDK
Symbol: SPANX
Synonyms: CT11.3; CTp11; SPANX-C; SPANX-E; SPANXE
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208858 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGACAAACAATCCAGTGCCGGCGGGGTGAAGAGGAGCGTCCCCTGTGAATCCAACGAGGTGAATGAGA
CGATGCCGGAGACCCCACTGGGGACTCAGACCCGCAACCTGCTCCTAAAAAATGAAACATCTGAGTC
CTCGACCATACTAGTGGTTCGCTACAGGAGGAACGTGAAAGAACATCTCCAGAGGAACCTGCTGAATGAC
CACGCCGAGAGAACAGAATCAACCCCTCCAATGGAGGAGGAGGAATTCATGGAAATAATGGTTGAAA
TACCTGCAAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208858 protein sequence
Red=Cloning site Green=Tags(s)

MDKQSSAGGVKRSVPCESNEVNETMPETPTGSDPQPAPKKMKTSESSTILVVRYYRRNVKRTSPEELLND
HARENRIPLQMEEEEFMEIMVEIPAK

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_022661
ORF Size: 291 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_022661.2, NP_073152.1</u>
RefSeq Size:	441 bp
RefSeq ORF:	294 bp
Locus ID:	64663
UniProt ID:	<u>Q9NY87</u>
Cytogenetics:	Xq27.2
MW:	11 kDa