

Product datasheet for RC208862A1V

Human CFHR3 (NM_021023) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CFHR3 (NM_021023) AAV Particle
Tag:	Myc-DDK
Symbol:	CFHR3
Synonyms:	CFHL3; DOWN16; FHR-3; FHR3; HLF4
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208862 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGTTGTTACTAATCAATGTCATTCTGACCTTGTGGGTTTCTGTGCTAATGGACAAGTGAACCTTGTC
 ATTTCCAGACATTAACATGGAGGTCTATTTTCATGAGAATATGCGTAGACCATACTTTCCAGTAGCTGT
 AGGAAAATATTACTCTATTACTGTGATGAACATTTTGAGACTCCTTCAGGAAGTTACTGGGATTACATT
 CATTGCACACAAAATGGGTGGTCACCAGCAGTACCATGTCTCAGAAAATGTTATTTTCCTTATTTGGAAA
 ATGGATATAATCAAAATTATGGAAGAAAGTTTGTACAGGGTAACCTACAGAAAGTTGCCATCCTGG
 CTACGGTCTTCCAAAAGCGCAGACCCAGTTACATGTACGGAGAAAGGCTGGTCTCCTACTCCCAGATGC
 ATCCGTGTCAGAACATGCTCAAAATCAGATATAGAAATTGAAAATGGATTCAATTCGGAATCTTCTCTA
 TTTATATTTTAAATAAAGAAATACAATATAAATGTAAACCAGGATATGCAACAGCAGATGGAAATCTTC
 AGGATCAATTACATGTTTGCAAAATGGATGGTCAGCACAACCAATTTGCATTAATCTTCAGAAAAGTGT
 GGGCTCCTCCACCTATTAGCAATGGTGATACCACTCTTTCTACTAAAAGTGTATGTGCCACAGTCAA
 GAGTCGAGTACCAATGCCAGCCCTACTATGAACCTCAGGGTTCTAATTATGTAAATGTAGTAATGGAGA
 GTGGTCGGAACCAAGATGCATACATCCATGTATAAATACTGAAGAAAACATGAATAAAAATAACATA
 AAGTTAAAGGAAGAAGTGACAGAAAATATTATGCAAAAACAGGGGATACCATTTGAATTTATGTGTAAT
 TGGGATATAATGCAAATACATCAATTCTATCATTTCAAGCAGTGTGTCGGAAGGGATAGTGAATACCC
 CAGATCGCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC208862 protein sequence Red=Cloning site Green=Tags(s) MLLLINVILTLWVSCANGQVKPCDFPDIKHGGLFHENMRRPYFPVAVGKYSSYYCDEHFETPSGSYWDYI HCTQNGWSPAVPCLRKCYFPYLENGYNQNYGRKFVQGNSTEVACHPGYGLPKAQTTVTCTEKGWSPTPRC IRVRTCSKSDIEIENGFISSSSIIYLNKEIQYKCKPGYATADGNSSGSITCLQNGWSAQFICINSSEKC GPPPIISNGDTSFLLKVYVPQSRVEYQCQPYELQGSNYVTCSNGEWSEPPRCIHPCIITEENMNKNNI KLKGRSDRKYYAKTGDTIEFMCKLGYNANTSILSFQAVCREGIVEYPRCE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_021023
ORF Size:	990 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_021023.4
RefSeq Size:	1667 bp
RefSeq ORF:	993 bp
Locus ID:	10878
UniProt ID:	Q02985
Cytogenetics:	1q31.3
MW:	37.3 kDa