

Product datasheet for **RC208872L1V**

Cathelicidin (CAMP) (NM_004345) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cathelicidin (CAMP) (NM_004345) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cathelicidin
Synonyms:	CAP-18; CAP18; CRAMP; FALL-39; FALL39; HSD26; LL37
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_004345
ORF Size:	510 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208872).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004345.3
RefSeq Size:	758 bp
RefSeq ORF:	513 bp
Locus ID:	820
UniProt ID:	P49913
Cytogenetics:	3p21.31
Domains:	Cathelicidins
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)

MW: 19.3 kDa

Gene Summary: This gene encodes a member of an antimicrobial peptide family, characterized by a highly conserved N-terminal signal peptide containing a cathelin domain and a structurally variable cationic antimicrobial peptide, which is produced by extracellular proteolysis from the C-terminus. In addition to its antibacterial, antifungal, and antiviral activities, the encoded protein functions in cell chemotaxis, immune mediator induction, and inflammatory response regulation. [provided by RefSeq, Sep 2014]