

Product datasheet for **RC208877A1V**

Human TTPA (NM_000370) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TTPA (NM_000370) AAV Particle
Tag:	Myc-DDK
Symbol:	TTPA
Synonyms:	alphaTTP; ATP; AVED; TTP1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208877 representing NM_000370 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCAGAGGCGCGATCCCAGCCCTCGGCGGGGCCGAGCTCAACGCGCTACCGGACCACTCTCCGTTGC
TGCAGCCGGGCTGGCGGCGCTGCGGCGCCGGGCCGGGAAGCTGGCGTCCCGCTCGCGCCGCTGCCGCT
CACCGACTCCTTCCTGCTGCGGTTCTGCGCGCCCGGATTCGATCTGGACCTGGCCTGGCGGTTACTA
AAAACTATTATAAGTGGAGAGCAGAATGTCCAGAAATAAGTGCAGATCTACACCCTAGAAGTATTATTG
GCCTCCTAAAGGCTGGCTACCATGGAGTCCTGAGATCCAGGGATCCCAGTGGCAGCAAGTCTTATTTA
CAGAATCGCACACTGGGACCCCAAAGTTTTTACAGCTTATGACGTATTTTCGAGTAAGTCTAATCACATCC
GAGCTTATTGTACAGGAGGTAGAACTCAGCGGAATGGAATCAAGGCTATCTTTGATCTGGAAGGTTGGC
AGTTTTCTCATGCTTTTCAAATCACTCCATCCGTAGCCAAGAAGATTGCTGCTGTACTTACGGATTCATT
TCCATTGAAAGTTCGTGGCATCCATTTGATAAATGAACCAGTAATTTTCATGCTGTCTTTCCATGATC
AAACCATTCTGACTGAAAAATTAAGGAACGGATTACATGCATGGGAACAACTACAAACAAAGCTTGC
TTCAGCATTTCCAGACATTCTCCTCTGGAATATGGTGGTGAAGAATTCTCATGGAGGACATTTGTCA
GGAATGGACAAATTTTATAATGAAGTCTGAAGATTATCTCAGCAGCATTTCTGAGAGCATTCATCA

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208877 representing NM_000370 Red=Cloning site Green=Tags(s) MAEARSQPSAGPQLNALPDHSPLLQPGLAALRRRAREAGVPLAPLPLTDSFLLRFLRARDFDLDLAWRL KNYYKWRAECPEISADLHPRSIIIGLLKAGYHGVLSRDPTGSKVLIYRIAHWDPKVFTAYDVFRVSLITS ELIVQEVETQRNGIKAIFDLEGWQFSHAFQITPSVAKKIAAVLTDSFPLKVRGIHLINPVIHAFVSMI KPFLTEKIKERIHMHGNNYQSLLQHFPDILPLEYGGEFSDMEDICQEWTFIMKSEDYLSSESISIQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_000370
ORF Size:	834 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_000370.2, NP_000361.1</u>
RefSeq Size:	1249 bp
RefSeq ORF:	837 bp
Locus ID:	7274
UniProt ID:	<u>P49638</u>
Cytogenetics:	8q12.3
MW:	31.6 kDa