

Product datasheet for **RC208885A1V**

Human EMG1 (NM_006331) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human EMG1 (NM_006331) AAV Particle
Tag:	Myc-DDK
Symbol:	EMG1
Synonyms:	C2F; Grcc2f; NEP1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC208885 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCGCGCCCAAGTGGATTCAAGCCTCGTGAACGAAGCGGTGGGAGCAGGCACAGGACTGGGATG
CTCTGCCACCCAAGCGGCCCCGACTAGGGGCAGGAAACAAGATCGGAGGCCGTAGGCTTATTGTGGTGCT
GGAAGGGGCCAGTCTGGAGACAGTCAAGGTAGGGAAGACATATGAGCTACTCAACTGTGACAAGCACAAG
TCTATATTGTTGAAGAATGGACGGGACCCTGGGAAGCGCGGCCAGATATCACCACCAAGAGTTTGCTGA
TGCTGATGGATAGTCCCCTGAACCGAGCTGGCTTGCTACAGGTTTATATCCATACACAGAAGAATGTTCT
GATTGAAGTGAATCCCCAGACCCGAATTCAGAACCTTTGACCGCTTTGTGGCCTCATGGTTCAACTT
TTACACAAGCTCAGTGTTTCGAGCAGCTGATGGCCCCAGAAGCTTTTGAAGTAATTAAGAATCCAGTAT
CAGATCACTTTCCAGTTGGATGTATGAAAGTTGGCACTTCTTTTCCATCCCGGTTGTCAGTGATGTGCG
TGAGCTGGTGCCAGCAGTATCCTATTGTTTTGTGGTAGGGGCCTTTGCCATGGCAAGGTCAGTGTC
GAGTATACAGAGAAGATGGTGTCCATCAGTAACACCCCTTTCTGCTGCCCTCACCTGTGCAAACTTA
CCACAGCCTTTGAGGAAGTATGGGGGTCATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC208885 protein sequence Red=Cloning site Green=Tags(s) MAAPSDGFKPRERSGGEQAQDWDALPPKRRLGAGNKIGRRRLIVLEGASLETVKVGKTYELLNCDKHK SILLKNGRDPGEARPDITHQSLLMLMDSPLNRAGLLQVYIHTQKNVLEIENPQTRIPRTFDRFCGLMVQL LHKLSVRAADGPQKLLKVIKNPVSDHFPVGCMKVGTSFSIPVSDVRELVPSSDPIVFVVGAFAGKVS EYTEKMVSISNYPLSAALTCAKLTAFEVWGVI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006331
ORF Size:	732 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006331.4</u>
RefSeq Size:	1072 bp
RefSeq ORF:	735 bp
Locus ID:	10436
UniProt ID:	<u>Q92979</u>
Cytogenetics:	12p13.31
MW:	26.7 kDa