

Product datasheet for **RC208885L4V**

EMG1 (NM_006331) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EMG1 (NM_006331) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EMG1
Synonyms:	C2F; Grcc2f; NEP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006331
ORF Size:	732 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208885).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006331.4
RefSeq Size:	1072 bp
RefSeq ORF:	735 bp
Locus ID:	10436
UniProt ID:	Q92979
Cytogenetics:	12p13.31
Domains:	Mra1
MW:	26.7 kDa



[View online »](#)

Gene Summary:

This gene encodes an essential, conserved eukaryotic protein that methylates pseudouridine in 18S rRNA. The related protein in yeast is a component of the small subunit processome and is essential for biogenesis of the ribosomal 40S subunit. A mutation in this gene has been associated with Bowen-Conradi syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]