

Product datasheet for **RC208890L4V**

SPATA19 (NM_174927) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SPATA19 (NM_174927) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SPATA19
Synonyms:	CT132; SPAS1; spergen1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_174927
ORF Size:	501 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208890).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_174927.1
RefSeq Size:	866 bp
RefSeq ORF:	504 bp
Locus ID:	219938
UniProt ID:	Q7Z5L4
Cytogenetics:	11q25
MW:	19 kDa
Gene Summary:	May have a role in spermiogenesis.[UniProtKB/Swiss-Prot Function]



[View online »](#)