

Product datasheet for RC208917L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MMP14 (NM_004995) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: MMP14 (NM_004995) Human Tagged ORF Clone Lentiviral Particle

Symbol: MMP14

Synonyms: MMP-14; MMP-X1; MT-MMP; MT-MMP 1; MT1-MMP; MT1MMP; MTMMP1; WNCHRS

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_004995 **ORF Size:** 1746 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC208917).

Sequence:
OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 004995.2

 RefSeq Size:
 3558 bp

 RefSeq ORF:
 1749 bp

 Locus ID:
 4323

 UniProt ID:
 P50281

 Cytogenetics:
 14q11.2

Domains: hemopexin, Peptidase_M10, ZnMc

Protein Families: Druggable Genome, Protease, Transmembrane







Protein Pathways: GnRH signaling pathway

MW: 65.89 kDa

Gene Summary: Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of

extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. However, the protein encoded by this gene is a member of the membrane-type MMP (MT-MMP) subfamily; each member of this subfamily contains a potential transmembrane domain suggesting that these proteins are expressed at the cell surface rather than secreted. This protein activates MMP2 protein, and this activity may be

involved in tumor invasion. [provided by RefSeq, Jul 2008]