

Product datasheet for **RC208929L4V**

FAM13B1 (FAM13B) (NM_016603) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM13B1 (FAM13B) (NM_016603) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM13B
Synonyms:	ARHGAP49; C5orf5; FAM13B1; KHCHP; N61
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016603
ORF Size:	2745 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208929).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016603.1
RefSeq Size:	5470 bp
RefSeq ORF:	2748 bp
Locus ID:	51306
UniProt ID:	Q9NYF5
Cytogenetics:	5q31.2
Domains:	RhoGAP
MW:	104.4 kDa



[View online »](#)