

Product datasheet for **RC208933L4V**

CTR9 (NM_014633) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CTR9 (NM_014633) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CTR9
Synonyms:	p150; p150TSP; SH2BP1; TSBP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014633
ORF Size:	3519 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208933).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014633.2
RefSeq Size:	4348 bp
RefSeq ORF:	3522 bp
Locus ID:	9646
UniProt ID:	Q6PD62
Cytogenetics:	11p15.4
Domains:	TPR
Protein Families:	Protease



[View online »](#)

MW: 133.3 kDa

Gene Summary: The protein encoded by this gene is a component of the PAF1 complex, which associates with RNA polymerase II and functions in transcriptional regulation and elongation. This complex also plays a role in the modification of histones. [provided by RefSeq, Oct 2016]