

Product datasheet for **RC208942L3V**

SREBP2 (SREBF2) (NM_004599) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SREBP2 (SREBF2) (NM_004599) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SREBP2
Synonyms:	bHLHd2; SREBP-2; SREBP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004599
ORF Size:	3423 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208942).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004599.2
RefSeq Size:	4325 bp
RefSeq ORF:	3426 bp
Locus ID:	6721
UniProt ID:	Q12772
Cytogenetics:	22q13.2
Domains:	HLH
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 123.5 kDa

Gene Summary: This gene encodes a member of the a ubiquitously expressed transcription factor that controls cholesterol homeostasis by regulating transcription of sterol-regulated genes. The encoded protein contains a basic helix-loop-helix-leucine zipper (bHLH-Zip) domain and binds the sterol regulatory element 1 motif. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]