

Product datasheet for **RC208944L1V**

SOX9 (NM_000346) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SOX9 (NM_000346) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SOX9 |
| Synonyms: | CMD1; CMPD1; SRA1; SRXX2; SRXY10 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_000346 |
| ORF Size: | 1527 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC208944). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_000346.2 |
| RefSeq Size: | 3935 bp |
| RefSeq ORF: | 1530 bp |
| Locus ID: | 6662 |
| UniProt ID: | P48436 |
| Cytogenetics: | 17q24.3 |
| Domains: | HMG |



[View online »](#)

| | |
|--------------------------|---|
| Protein Families: | Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Transcription Factors |
| MW: | 56 kDa |
| Gene Summary: | The protein encoded by this gene recognizes the sequence CCTTGAG along with other members of the HMG-box class DNA-binding proteins. It acts during chondrocyte differentiation and, with steroidogenic factor 1, regulates transcription of the anti-Muellerian hormone (AMH) gene. Deficiencies lead to the skeletal malformation syndrome campomelic dysplasia, frequently with sex reversal. [provided by RefSeq, Jul 2008] |