

Product datasheet for **RC208952L3V**

IKIP (IKBIP) (NM_153687) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IKIP (IKBIP) (NM_153687) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IKIP
Synonyms:	IKIP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153687
ORF Size:	1131 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208952).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153687.2
RefSeq Size:	3175 bp
RefSeq ORF:	1134 bp
Locus ID:	121457
UniProt ID:	Q70UQ0
Cytogenetics:	12q23.1
MW:	43.1 kDa
Gene Summary:	Target of p53/TP53 with pro-apoptotic function.[UniProtKB/Swiss-Prot Function]


[View online »](#)