

## Product datasheet for RC208957L4V

### SKAR (POLDIP3) (NM\_178136) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SKAR (POLDIP3) (NM_178136) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SKAR
Synonyms:	PDIP3; PDIP46; SKAR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178136
ORF Size:	1176 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208957).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178136.1</a>
RefSeq Size:	3348 bp
RefSeq ORF:	1179 bp
Locus ID:	84271
UniProt ID:	<a href="#">Q9BY77</a>
Cytogenetics:	22q13.2
MW:	42.7 kDa



[View online »](#)

**Gene Summary:**

This gene encodes an RRM (RNA recognition motif)-containing protein that participates in the regulation of translation by recruiting ribosomal protein S6 kinase beta-1 to mRNAs. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]