

Product datasheet for **RC208959A1V**

Human **LRRC57 (NM_153260)** AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human LRRC57 (NM_153260) AAV Particle
Tag: Myc-DDK
Symbol: LRRC57
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC208959 representing NM_153260
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGAAACAGTGCCTCCGCGCTCATGTGAAACGGCGCAGAAAACGGTGTCTTTCAGCTTAAGGACC
GAGGGCTGACCGAGTTCCTCCGAGACTGCAGAAGCTGACGAGCAATCTCAGGACCATCGACTTGCCAA
CAACAAGATCGAAAGCCTACCGCCTTTGCTGATAGGAAAGTTCACCTGCTGAAGAGCCTCTCCCTGAAC
AACAACTGACTGTTCTGCCTGATGAGATATGCAATCTGAAAAACTAGAGACGCTAAGCCTAAACA
ACAATCACCTTAGAGAGCTGCCGTCTACCTTTGGGCAACTCTCTGCCCTCAAGACCTGAGCCTCTCTGG
GAACCAACTGGGAGCATTACCTCCCAACTTTGTAGCCTACGGCACCTGGATGTGATGGATCTCTCTAAG
AACCAGATTGCAAGTATACCTGACTCAGTGGGAGAGCTGCAAGTCATCGAACTCAACCTCAACGAGAATC
AGATATCTCAGATCTCAGTGAAGATATCTTGTGTCCACGCCTTAAAAATCTTCGCTGGAAGAGAATTG
TCTTGAGCTCAGCATGCTTCCCGAGAGCATCTCAGTATTCCAGATCTGTCTGCTGTGGAAGGC
AATCTTTTGAATAAAGAACTTCGAGAACTGGAAGGCTATGATAAGTACATGGAGAGGTTACAGCCA
CCAAGAAGAAGTTTGCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC208959 representing NM_153260
Red=Cloning site Green=Tags(s)

MGNSALRAHVETAQKTGVFQLKDRGLTEFPADLQKLTSNLRITDLSNNKIESLPPLLIGKFTLLKSLSLN
NNKLTVPDEICNLKKLETLSLNNHLRELPTFGQLSALKTSLSGNQLGALPPQLCSLRHLDVMDLSK
NQIRIPDSVGEIQVIELNLNQNISQISVKISCCPRLKILRLEENCLELSMLPQSILSDSQICLLAVEG
NLFEIKKLRELEGYDKYMERFTATKKKFA

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_153260
ORF Size:	717 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_153260.1
RefSeq Size:	2647 bp
RefSeq ORF:	720 bp
Locus ID:	255252
UniProt ID:	Q8N9N7
Cytogenetics:	15q15.2
MW:	26.6 kDa