

## Product datasheet for **RC208968L4V**

### PSKH1 (NM\_006742) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PSKH1 (NM_006742) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PSKH1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006742
ORF Size:	1272 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208968).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006742.1</a>
RefSeq Size:	3505 bp
RefSeq ORF:	1275 bp
Locus ID:	5681
UniProt ID:	<a href="#">P11801</a>
Cytogenetics:	16q22.1
Protein Families:	Druggable Genome, Protein Kinase
MW:	48 kDa



[View online »](#)

**Gene Summary:**

May be a SFC-associated serine kinase (splicing factor compartment-associated serine kinase) with a role in intranuclear SR protein (non-snRNP splicing factors containing a serine/arginine-rich domain) trafficking and pre-mRNA processing.[UniProtKB/Swiss-Prot Function]