

## Product datasheet for **RC208974L2V**

### USP38 (NM\_032557) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	USP38 (NM_032557) Human Tagged ORF Clone Lentiviral Particle
Symbol:	USP38
Synonyms:	HP43.8KD
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_032557
ORF Size:	3126 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208974).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032557.4</a>
RefSeq Size:	7081 bp
RefSeq ORF:	3129 bp
Locus ID:	84640
UniProt ID:	<a href="#">Q8NB14</a>
Cytogenetics:	4q31.21
Domains:	UCH
Protein Families:	Druggable Genome, Protease



[View online »](#)

**MW:** 116.5 kDa

**Gene Summary:** Deubiquitinating enzyme exhibiting a preference towards 'Lys-63'-linked ubiquitin chains.  
[UniProtKB/Swiss-Prot Function]