

Product datasheet for **RC208991L1V**

KREMEN1 (NM_001039571) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KREMEN1 (NM_001039571) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KREMEN1
Synonyms:	FLJ31863; KREMEM1; KREMEN; kringle-coding gene marking the eye and the nose; kringle-containing transmembrane protein 1; kringle containing transmembrane protein 1; KRM1; OTTHUMP00000028977
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001039571
ORF Size:	1374 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC208991).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001039571.1 , NP_001034660.1
RefSeq Size:	6115 bp
RefSeq ORF:	1376 bp
Locus ID:	83999
Cytogenetics:	22q12.1
Protein Families:	Druggable Genome, Transmembrane
MW:	50.18 kDa



[View online »](#)

Gene Summary:

This gene encodes a high-affinity dickkopf homolog 1 (DKK1) transmembrane receptor that functionally cooperates with DKK1 to block wingless (WNT)/beta-catenin signaling. The encoded protein is a component of a membrane complex that modulates canonical WNT signaling through lipoprotein receptor-related protein 6 (LRP6). It contains extracellular kringle, WSC, and CUB domains. Alternatively spliced transcript variants encoding distinct isoforms have been observed for this gene. [provided by RefSeq, Jul 2008]