

Product datasheet for **RC209026L4V**

CTDSP2 (NM_005730) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CTDSP2 (NM_005730) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CTDSP2
Synonyms:	OS4; PSR2; SCP2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005730
ORF Size:	813 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209026).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005730.3
RefSeq Size:	5038 bp
RefSeq ORF:	816 bp
Locus ID:	10106
UniProt ID:	O14595
Cytogenetics:	12q14.1
Domains:	CPDc
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)

MW: 30.7 kDa

Gene Summary: Preferentially catalyzes the dephosphorylation of 'Ser-5' within the tandem 7 residue repeats in the C-terminal domain (CTD) of the largest RNA polymerase II subunit POLR2A. Negatively regulates RNA polymerase II transcription, possibly by controlling the transition from initiation/capping to processive transcript elongation. Recruited by REST to neuronal genes that contain RE-1 elements, leading to neuronal gene silencing in non-neuronal cells. May contribute to the development of sarcomas.[UniProtKB/Swiss-Prot Function]