

Product datasheet for **RC209070L4V**

Phospholipase D1 (PLD1) (NM_002662) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Phospholipase D1 (PLD1) (NM_002662) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PLD1
Synonyms:	CVDD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002662
ORF Size:	3222 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209070).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002662.2
RefSeq Size:	6025 bp
RefSeq ORF:	3225 bp
Locus ID:	5337
UniProt ID:	Q13393
Cytogenetics:	3q26.31
Protein Pathways:	Endocytosis, Ether lipid metabolism, Fc gamma R-mediated phagocytosis, Glycerophospholipid metabolism, GnRH signaling pathway, Metabolic pathways, Pancreatic cancer, Pathways in cancer



[View online »](#)

MW: 124.2 kDa

Gene Summary: This gene encodes a phosphatidylcholine-specific phospholipase which catalyzes the hydrolysis of phosphatidylcholine in order to yield phosphatidic acid and choline. The enzyme may play a role in signal transduction and subcellular trafficking. Alternative splicing results in multiple transcript variants with both catalytic and regulatory properties. [provided by RefSeq, Sep 2011]