

Product datasheet for **RC209072L3V**

Chymotrypsin like protease (CTRL) (NM_001907) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Chymotrypsin like protease (CTRL) (NM_001907) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Chymotrypsin like protease
Synonyms:	CTRL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001907
ORF Size:	792 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209072).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001907.1
RefSeq Size:	1197 bp
RefSeq ORF:	795 bp
Locus ID:	1506
UniProt ID:	P40313
Cytogenetics:	16q22.1
Protein Families:	Druggable Genome, Protease



[View online »](#)

MW: 28 kDa

Gene Summary: This gene encodes a serine-type endopeptidase with chymotrypsin- and elastase-2-like activities. The gene encoding this zymogen is expressed specifically in the pancreas and likely functions as a digestive enzyme. [provided by RefSeq, Sep 2016]