

Product datasheet for **RC209075L1V**

TUBB2B (NM_178012) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | TUBB2B (NM_178012) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | TUBB2B |
| Synonyms: | bA506K6.1; CDCBM7; PMGYSA |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_178012 |
| ORF Size: | 1335 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC209075). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_178012.3 |
| RefSeq Size: | 2019 bp |
| RefSeq ORF: | 1338 bp |
| Locus ID: | 347733 |
| UniProt ID: | Q9BVA1 |
| Cytogenetics: | 6p25.2 |
| Protein Families: | Druggable Genome |
| Protein Pathways: | Gap junction, Pathogenic Escherichia coli infection |



[View online »](#)

MW: 50 kDa

Gene Summary: The protein encoded by this gene is a beta isoform of tubulin, which binds GTP and is a major component of microtubules. This gene is highly similar to TUBB2A and TUBB2C. Defects in this gene are a cause of asymmetric polymicrogyria. [provided by RefSeq, Mar 2010]