

Product datasheet for **RC209098A1V**

Human HBXIP (LAMTOR5) (NM_006402) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human HBXIP (LAMTOR5) (NM_006402) AAV Particle
Tag: Myc-DDK
Symbol: HBXIP
Synonyms: HBXIP; XIP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC209098 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGCCAGGTGCAGGTCACCTCGACGGTCACCGCGCGGGGAGCCCAAGCCTTCGTCAGGCTCTGTGCG
 ACGGAAGCGCAGTGATGTTTTCCAGTAAAGAACGCGGACGTTGCACCGTGATCAATTTTGTCCCTTTGGA
 GGCGCCGTTACGGTCCACGCCCGCTCGCGTCAAGTACTGAGGCCTGTGGTGGAGAAGGACGTGCCGTG
 CCGCTGGGTTCTGAGCCGGAGTGGTGGTGGTGGGATGGAGGCGACCTTGGAGCAGCACTTGAAGACA
 CAATGAAGAATCCCTCCATTGTTGGAGTCTGTGCACAGATTACAAGGACTTAATCTGGGTTGCCGCGG
 GACCCTGTGAGATGAGCATGCTGGAGTGATATCTGTTCTAGCCAGCAAGCAGCTAAGCTAACCTCTGAC
 CCCACTGATATTCCTGTGGTGTGTCTAGAATCAGATAATGGGAACATTATGATCCAGAAACACGATGGCA
 TCACGGTGGCAGTGCACAAAATGGCCTCT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC209098 protein sequence
 Red=Cloning site Green=Tags(s)

MEPGAGHLDGHRAGSPSLRQALCDGSAYMFSSKERGRCTVINFPLEAPLRSTPRSRQVTEACGGEGRAV
 PLGSEPEWSVGMEATLEQHLEDTMKNPSIVGLCTDSQGLNLGCRGTLSDHAGVISVLAQQAALTS
 PTDIPVVCLESDNGNIMIQLKHDGITVAVHKMAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_006402
ORF Size:	519 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006402.2</u>
RefSeq Size:	855 bp
RefSeq ORF:	522 bp
Locus ID:	10542
UniProt ID:	<u>O43504</u>
Cytogenetics:	1p13.3
MW:	18.2 kDa