

Product datasheet for **RC209137L1V**

GRK5 (NM_005308) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GRK5 (NM_005308) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GRK5
Synonyms:	FP2025; GPRK5
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005308
ORF Size:	1770 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209137).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005308.2
RefSeq Size:	2575 bp
RefSeq ORF:	1773 bp
Locus ID:	2869
UniProt ID:	P34947
Cytogenetics:	10q26.11
Domains:	RGS, pkinase, S_TK_X, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase


[View online »](#)

Protein Pathways: Chemokine signaling pathway, Endocytosis

MW: 67.8 kDa

Gene Summary: This gene encodes a member of the guanine nucleotide-binding protein (G protein)-coupled receptor kinase subfamily of the Ser/Thr protein kinase family. The protein phosphorylates the activated forms of G protein-coupled receptors thus initiating their deactivation. It has also been shown to play a role in regulating the motility of polymorphonuclear leukocytes (PMNs). [provided by RefSeq, Jul 2008]