

Product datasheet for **RC209163L3V**

AGO1 (NM_012199) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AGO1 (NM_012199) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AGO1
Synonyms:	EIF2C; EIF2C1; GERP95; hAgo1; Q99
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012199
ORF Size:	2571 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209163).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012199.2
RefSeq Size:	7478 bp
RefSeq ORF:	2574 bp
Locus ID:	26523
UniProt ID:	Q9UL18
Cytogenetics:	1p34.3
Protein Families:	Druggable Genome
MW:	97 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the argonaute family of proteins, which associate with small RNAs and have important roles in RNA interference (RNAi) and RNA silencing. This protein binds to microRNAs (miRNAs) or small interfering RNAs (siRNAs) and represses translation of mRNAs that are complementary to them. It is also involved in transcriptional gene silencing (TGS) of promoter regions that are complementary to bound short antigenic RNAs (agRNAs), as well as in the degradation of miRNA-bound mRNA targets. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. A recent study showed this gene to be an authentic stop codon readthrough target, and that its mRNA could give rise to an additional C-terminally extended isoform by use of an alternative in-frame translation termination codon. [provided by RefSeq, Nov 2015]