

## Product datasheet for **RC209168L2V**

### PGM2L1 (NM\_173582) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PGM2L1 (NM_173582) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PGM2L1
Synonyms:	BM32A; PMMLP
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_173582
ORF Size:	1866 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209168).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_173582.3</a>
RefSeq Size:	8519 bp
RefSeq ORF:	1869 bp
Locus ID:	283209
UniProt ID:	<a href="#">Q6PCE3</a>
Cytogenetics:	11q13.4
Protein Pathways:	Starch and sucrose metabolism
MW:	70.4 kDa



[View online »](#)

**Gene Summary:**

Glucose 1,6-bisphosphate synthase using 1,3-bisphosphoglycerate as a phosphate donor and a series of 1-phosphate sugars as acceptors, including glucose 1-phosphate, mannose 1-phosphate, ribose 1-phosphate and deoxyribose 1-phosphate. 5 or 6-phosphosugars are bad substrates, with the exception of glucose 6-phosphate. Also synthesizes ribose 1,5-bisphosphate. Has only low phosphopentomutase and phosphoglucomutase activities. [UniProtKB/Swiss-Prot Function]