

Product datasheet for **RC209171L4V**

Centaurin beta 2 (ACAP2) (NM_012287) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Centaurin beta 2 (ACAP2) (NM_012287) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Centaurin beta 2
Synonyms:	CENTB2; CNT-B2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_012287
ORF Size:	2334 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209171).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012287.3 , NP_036419.2
RefSeq Size:	7187 bp
RefSeq ORF:	2337 bp
Locus ID:	23527
UniProt ID:	Q15057
Cytogenetics:	3q29
Domains:	ArfGap, PH, ANK
Protein Pathways:	Endocytosis


[View online »](#)

MW: 88.1 kDa

Gene Summary: GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6).[UniProtKB/Swiss-Prot Function]