

Product datasheet for **RC209177L3V**

ST6GALNAC3 (NM_152996) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ST6GALNAC3 (NM_152996) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ST6GALNAC3
Synonyms:	PRO7177; SIAT7C; ST6GALNACIII; STY
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_152996
ORF Size:	915 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209177).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152996.1
RefSeq Size:	3259 bp
RefSeq ORF:	918 bp
Locus ID:	256435
UniProt ID:	Q8NDV1
Cytogenetics:	1p31.1
Domains:	Glyco_transf_29
Protein Families:	Transmembrane



[View online »](#)

Protein Pathways: Glycosphingolipid biosynthesis - ganglio series, Metabolic pathways

MW: 35.4 kDa

Gene Summary: ST6GALNAC3 belongs to a family of sialyltransferases that transfer sialic acids from CMP-sialic acid to terminal positions of carbohydrate groups in glycoproteins and glycolipids (Lee et al., 1999 [PubMed 10207017]).[supplied by OMIM, Mar 2008]