

Product datasheet for **RC209219A1V**

Human TRIB1 (NM_025195) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TRIB1 (NM_025195) AAV Particle
Tag:	Myc-DDK
Symbol:	TRIB1
Synonyms:	C8FW; GIG-2; GIG2; SKIP1; TRB-1; TRB1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC209219 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCGGGTCCGGTCCGGTGCCTCTGCCATGAGCGGCGCCTCGCAGCCCCGCGGCCCGGCCCTGCTCTTCC
CAGCCACCCGAGGCGTCCCGGCCAAACGCTGCTGGACGCCGACGACGCGCGGCTGTGGCGGCCAAGTG
CCCGCGCTCTCCGAGTGCTCCAGCCCCCGGACTACCTCAGCCCCCGGCTCGCCCTGCAGTCCGCAG
CCCCCGCTGCCGCTCCGGGGCCGGCGGAGGCTCCGGGAGCGCGCCGGGGCCAGCCGCATCGCCGACT
ACCTGCTGCTGCCCTAGCCGAGCGCGAGCATGTGTCCGGGCGCTGTGCATCCACACTGGACGCGAGCT
GCGCTGCAAGGTGTTTCCATTAAACACTACCAGGACAAAATCAGGCCTTACATCCAGCTGCCATCGCAC
AGCAACATTACTGGCATTGTGGAAGTGATCCTTGGGAAACCAAGGCCTATGTCTTTGAGAAGGACT
TTGGGGACATGCACTCCTATGTGCGAAGCCGGAAGAGGCTGCGGGAAGAGGAAGCCGCCGCTCTTCAA
GCAGATTGTCTCCGCGTCCGCCACTGCCACCAGTCAGCCATCGTGTGGGGACCTGAAGCTTAGGAAG
TTCGTCTTCCACGGAGGAGAGAACCAGCTTAGACTAGAAAGTCTAGAAGACACACATAATGAAGG
GGGAAGATGATGCTTTGTCAGACAAACATGGCTGCCAGCCTACGTGAGCCCTGAGATCTCAACACCAC
TGGGACCTACTCCGAAAGGCTGCGGACGTTTGGAGCCTGGGGTGATGCTCTACACCCTCTGGTTGGA
CGATACCCCTTCCATGACTCAGACCCAGTGCCCTTTTCTCAAATTCGGCGTGACAGTTCTGCATTC
CTGAGCACATTTCCCAAGGAGGTCCTCATTGCGAGCCTTTGAGACGGGAGCCCTCCGAGAGACT
CACTGCCCCGAGATCCTACTGCACCCCTGTTTGTAGTCCGTCTTGAACCCGGGTACATCGACTCAGAA
ATAGGAACTTCAGACCAGATTGTTCCAGAGTACCAGGAGGACAGTGACATTAGTTCCTTCTTCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>RC209219 protein sequence Red=Cloning site Green=Tags(s) MRVGPVRSAMSGASQPRGPALLFPATRGVPAKRLLDADDAVAACCPRLSECSSPPDYLSPPGSPCSPQ PPPAAPGAGGSGSAPGPSRIADYLLPLAEREHVSRALCIHTGRELRCKVFPIKHYQDKIRPYIQLPSH SNITGIVEVILGETKAYVFFEKDFGDMHSYVRSRKRLREEEAARLFKQIVSAVAHCHQSAIVLGDLKLRK FVFSTEERTQLRLESLEDTHIMKGEDDALSDKHGCPAYVSPEILNTTGTYSGKAADVWSLGVMLYTLLVG RYPFHSDPSALFSKIRRGQFCIPEHISPKARCLIRSLRREPSERLTAPEILLHPWFESVLEPGYIDSE IGTSDQIVPEYQEDSDISSFFC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_025195
ORF Size:	1116 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_025195.2
RefSeq Size:	3649 bp
RefSeq ORF:	1119 bp
Locus ID:	10221
UniProt ID:	Q96RU8
Cytogenetics:	8q24.13
MW:	41 kDa