

Product datasheet for **RC209222L1V**

ENPP1 (NM_006208) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ENPP1 (NM_006208) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ENPP1
Synonyms:	ARHR2; COLED; M6S1; NPP1; NPPS; PC-1; PCA1; PDNP1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_006208
ORF Size:	2619 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209222).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006208.1 , NP_006199.1
RefSeq Size:	7442 bp
RefSeq ORF:	2778 bp
Locus ID:	5167
UniProt ID:	P22413
Cytogenetics:	6q23.2
Domains:	SO, Endonuclease, Phosphodiester
Protein Families:	Druggable Genome, Transmembrane




[View online »](#)

Protein Pathways: Metabolic pathways, Nicotinate and nicotinamide metabolism, Pantothenate and CoA biosynthesis, Purine metabolism, Riboflavin metabolism, Starch and sucrose metabolism

MW: 99.9 kDa

Gene Summary: This gene is a member of the ecto-nucleotide pyrophosphatase/phosphodiesterase (ENPP) family. The encoded protein is a type II transmembrane glycoprotein comprising two identical disulfide-bonded subunits. This protein has broad specificity and cleaves a variety of substrates, including phosphodiester bonds of nucleotides and nucleotide sugars and pyrophosphate bonds of nucleotides and nucleotide sugars. This protein may function to hydrolyze nucleoside 5' triphosphates to their corresponding monophosphates and may also hydrolyze diadenosine polyphosphates. Mutations in this gene have been associated with 'idiopathic' infantile arterial calcification, ossification of the posterior longitudinal ligament of the spine (OPLL), and insulin resistance. [provided by RefSeq, Jul 2008]

Product images:



[RC209222L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC209222L1V particle to overexpress human ENPP1-Myc-DDK fusion protein.