

Product datasheet for **RC209238L4V**

AGTPBP1 (NM_015239) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	AGTPBP1 (NM_015239) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AGTPBP1
Synonyms:	CCP1; CONDCA; NNA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015239
ORF Size:	3558 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209238).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015239.1
RefSeq Size:	4436 bp
RefSeq ORF:	3561 bp
Locus ID:	23287
UniProt ID:	Q9UPW5
Cytogenetics:	9q21.33
Domains:	Zn_carbOpept
Protein Families:	Protease



[View online »](#)

MW: 133.9 kDa

Gene Summary: NNA1 is a zinc carboxypeptidase that contains nuclear localization signals and an ATP/GTP-binding motif that was initially cloned from regenerating spinal cord neurons of the mouse. [supplied by OMIM, Jul 2002]