

Product datasheet for **RC209240A1V**

Human KIAA0247 (SUSD6) (NM_014734) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human KIAA0247 (SUSD6) (NM_014734) AAV Particle
Tag:	Myc-DDK
Symbol:	KIAA0247
Synonyms:	DRAGO; KIAA0247
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC209240 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGCCATGGCAGGATAGCACCAAAGAGCACCTCAGTGTTTGCCGTGGCCTCCGTGGGACATGGAGTGT
TCCTTCCGCTAGTGATCCTTTGCACCCTGCTTGGAGACGGACTTGCTTCCGTGTGCCCTACCACCGGA
GCCAGAGAATGGTGGCTACATCTGCCACCCCGGCCCTGCAGAGACCCCTGACAGCAGGCAGTGTCATC
GAATACCTGTGTGCTGAAGGCTACATGTTGAAGGGCGATTACAAATACCTGACGTGTAAGAATGGCGAGT
GGAAACCAGCCATGGAGATTAGCTGCCGTCTCAACGAGGATAAAGACACCCACACATCACTTGGGGTCCC
CACGCTGTCTATAGTGGCTTCTACTGCCAGCTCCGTGGCGCTCATTCTCCTCCTCGTGGTGCTTTGTG
CTGCTGCAGCCAAAGCTGAAGTCTTTCCATCATAGCAGGCGTGACCAGGGGTATCTGGGACCAGGTCT
CCATCATGGTGGATGGAGTCCAGGTTGCACTACCATCATACGAGGAGGCTGTATATGGCAGTTCTGGTCA
CTGTGTGCCACCTGCTGACCCAGAGTACAGATTGTGCTGTCAGAAGGGTCTGGGCCAGTGGGAGGAGC
GTGCCAAGGGAGCAACAGCTGCCGACCAAGGGCCTGCTCCTCTGCAGGTGGAGAAGATGAGGCCCCAG
GCCAGTCTGGACTATGTGAAGCCTGGGGCTCTCGGGCCTCAGAGACTGTGATGGTGCATCAGGCAACCAC
CTCTTCTGGGTGGCCGGCTCAGGGAACCGCAACTGGCACACAAAGAACTGCAGATTCAGAGAACAGT
GACATACAAAGCCTTTTATCCCTCACGTGAGGAGTACACAGATGATATCCACTGTTGAAAGAAGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC209240 protein sequence Red=Cloning site Green=Tags(s) MCHGRIAPKSTSVFAVASVGHGVFLPLVILCTLLGDGLASVCPLPPEPENGGYICHPRPCRDPLTAGSVI EYLCAEGYMLKGDYKYLTCNGEWKPAMEISCRLNEDKDTHTSLGVPTLSIVASTASSVALILLVVLV LLQPKLKSFFHHSRRDQGVSGDQVSIMVDGVQVALPSYEEAVYGSSGHCVPPADPRVQIVLSESGPSGRS VPREQLPDQGACSSAGGEDEAPGQSGLCEAWGSRASETVMVHQATTSSWVAGSGNRQLAHKETADSENS DIQSLLSLTSEEYTDIPLLKEA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_014734
ORF Size:	909 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_014734.2
RefSeq Size:	5393 bp
RefSeq ORF:	912 bp
Locus ID:	9766
UniProt ID:	Q92537
Cytogenetics:	14q24.1
MW:	32.1 kDa