

## Product datasheet for **RC209258L2V**

### Niemann Pick C1 (NPC1) (NM\_000271) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Niemann Pick C1 (NPC1) (NM_000271) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NPC1
Synonyms:	NPC; POGZ; SLC65A1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_000271
ORF Size:	3834 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209258).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000271.4</a>
RefSeq Size:	4827 bp
RefSeq ORF:	3837 bp
Locus ID:	4864
UniProt ID:	<a href="#">O15118</a>
Cytogenetics:	18q11.2
Domains:	Patched
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**Protein Pathways:** Lysosome

**MW:** 142.1 kDa

**Gene Summary:** This gene encodes a large protein that resides in the limiting membrane of endosomes and lysosomes and mediates intracellular cholesterol trafficking via binding of cholesterol to its N-terminal domain. It is predicted to have a cytoplasmic C-terminus, 13 transmembrane domains, and 3 large loops in the lumen of the endosome - the last loop being at the N-terminus. This protein transports low-density lipoproteins to late endosomal/lysosomal compartments where they are hydrolyzed and released as free cholesterol. Defects in this gene cause Niemann-Pick type C disease, a rare autosomal recessive neurodegenerative disorder characterized by over accumulation of cholesterol and glycosphingolipids in late endosomal/lysosomal compartments.[provided by RefSeq, Aug 2009]