

Product datasheet for **RC209274L3V**

HAPLN1 (NM_001884) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HAPLN1 (NM_001884) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HAPLN1
Synonyms:	CRT1; CRTL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001884
ORF Size:	1062 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209274).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001884.2
RefSeq Size:	4678 bp
RefSeq ORF:	1065 bp
Locus ID:	1404
UniProt ID:	P10915
Cytogenetics:	5q14.3
Domains:	Xlink, ig, IGv, IG
MW:	40.1 kDa



[View online »](#)

Gene Summary:

Stabilizes the aggregates of proteoglycan monomers with hyaluronic acid in the extracellular cartilage matrix.[UniProtKB/Swiss-Prot Function]