

Product datasheet for **RC209280A1V**

Human TMEM88 (NM_203411) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: TMEM88

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC209280 representing NM_203411
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGGATGTCCCGGGGCACAGCGAGCGGTTCTGGTGACGGCCAGAGCCCGGGACCCCTGGACT
GTTGGGCTGCGCTGTTCTTGTAAACAGCCAGAATCTGCTGGTGGCTGCCTTCAATCTTCTCCTGCTGGT
GCTGGTGCTAGGGACCATCTTGCTACCCGCTGTCACCATGCTGGGCTTCGGCTTCTCTGCCACTCTCAG
TTCCTGCGCTCCAGGCACCCCTTGACCGCGCACCTGCGGGACCCGGTTTCACGGCCCTACTGGTCA
CCGGATTCTGCTCCTCGTGCCGCTGCTCGTGCTTGCTCTGGCCAGCTACCGCCGCTCTGCCTGCGCCT
CCGCTAGCCGATTGCCTCGTGCCCTACAGCCGAGCCCTTTATCGGCGTCGGCGCGCCCGCAGCCGCGG
CAATCCGGGCTCACCAGGTCACAGGCCGTTCCACATCAGGAAAGGTCTGGGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC209280 representing NM_203411
Red=Cloning site Green=Tags(s)

MADVPGAQRAVPGDGPEDPLDCWACAVLVTAQNLLVAAFNLLLLVLVLTILLPAVMTLGFGLCHSQ
FLRSQAPPCTAHLRDPGFTALLVTGFLLLVPLLVLALASYRRLCLRLADCLVPYSRALYRRRRAPQPR
QIRASPGSQAVPTSGKVWV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_203411



ORF Size:	477 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_203411.1
RefSeq Size:	887 bp
RefSeq ORF:	480 bp
Locus ID:	92162
UniProt ID:	Q6PEY1
Cytogenetics:	17p13.1
MW:	17.1 kDa