

## Product datasheet for **RC209295L3V**

### **BCL6B (NM\_181844) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BCL6B (NM_181844) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BCL6B
Synonyms:	BAZF; ZBTB28; ZNF62
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_181844
ORF Size:	1440 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209295).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_181844.2</a>
RefSeq Size:	3544 bp
RefSeq ORF:	1440 bp
Locus ID:	255877
UniProt ID:	<a href="#">Q8N143</a>
Cytogenetics:	17p13.1
MW:	51.6 kDa



[View online »](#)

**Gene Summary:**

Acts as a sequence-specific transcriptional repressor in association with BCL6. May function in a narrow stage or be related to some events in the early B-cell development.  
[UniProtKB/Swiss-Prot Function]