

Product datasheet for RC209296

C1orf122 (NM_198446) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C1orf122 (NM_198446) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C1orf122
Synonyms: ALAESM
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC209296 representing NM_198446
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAATGGGGCCCGGGCTCAGACTGGTCACGGGGGAGGCTGCCGGCGTGGACCGGGAAAGCGGGGC
 TGGGGCTCGGCGGGAGGCCACCCACAGCCGCCCGGGAGGAGCGCCAGCAGCTGCTGGACGCGGT
 GGAGCAGCGGCAGCGGCAGCTCCTGGACACCATCGCAGCCTGCGAGGAGATGTTACGGCAGCTGGGCCGC
 CGGCGCCCGAGCCGGCTGGTGGCGGGAACGTCTCAGCCAAACCTGGAGCGCCCCCAGCCGGCTGCT
 CGCCAGAGGCGGCTTCCAAAGGATGTGGCGATGGAGCTGCGGAGCCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC209296 representing NM_198446
 Red=Cloning site Green=Tags(s)
 MEWGPGSDWSRGEAAGVDRGKAGLGLGGRPPPQPPREERAQQLLDAVEQRQRQLLDLIAACEMLRQLGR
 RRPEPAGGGNVSAPGAPPQPAVSARGGFPKDAGDGAEP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_198446

ORF Size: 330 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


RefSeq: [NM_198446.3](#)

RefSeq Size: 1359 bp

RefSeq ORF: 333 bp

Locus ID: 127687
 UniProt ID: [Q6ZSJ8](#)
 Cytogenetics: 1p34.3
 MW: 11.5 kDa

Product images:



Circular map for RC209296