

Product datasheet for **RC209315A1V**

Human ZDHHC24 (NM_207340) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: ZDHHC24

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC209315 representing NM_207340
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGGGCAGCCCTGGGCGCTGGGAGCACGGACGGGGCGCCGCGCAGCTGCCTCTCGTGCTCACCGCGC
TGTGGGCCGCGGCCGTGGGCCTGGAGCTGGCTTACGTGCTGGTGCTCGGTCCCGGGCCGCCCGCTGGG
ACCCCTGGCCCGGCCCTGCAGCTGGCGCTGGCCGCCTTCCAGCTGCTCAACCTGCTGGGCAACGTGGG
CTCTTCTGCGCTCGGATCCCAGCATCCGTGGCGTGATGCTGGCCGGCCGCGTCTGGGCCAGGGCTGGG
CTTACTGCTACCAATGCCAAAGCCAGGTGCCGCCACGCAGCGGACACTGCTCTGCCTGCCGCTCTGCAT
CCTGCGTCGGGACCACTGCCGCTGCTGGGCGCTGCGTGGGCTTCGGCAACTACCGGCCCTTCCTG
TGCTGCTGCTTCATGCCGCCGGGCTCCTGCTCCACGTCTCTGTGCTGCTGGGCCCTGCACTGTCGGCCC
TGCTGCGAGCCACAGCCCCCTCCACATGGCTGCCCTCCTCCTGCTCCCTGGCTCATGTTGCTCACAGG
CAGAGTGTCTCTGGCACAGTTTGCCTTGGCCTTCGTGACGGACACGTGCGTGGCGGGTGCGCTGCTGTG
GGGGCTGGGCTGCTCTTCATGGGATGCTGCTGCTGCGGGGCCAGACCACATGGGAGTGGGCTCGGGGCC
AGCACTCCTATGACCTGGGTCCTGCCACAACCTGCAGGCAGCCCTGGGGCCCCGCTGGGCCCTCGTCTG
GCTCTGGCCCTTCTGGCCTCCCCATTGCCTGGGGATGGGATCACCTTCAGACCACAGCAGATGTGGGA
CACACAGCCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>RC209315 representing NM_207340 Red=Cloning site Green=Tags(s) MGQPWAAAGSTDGAPQLPLVLTALWAAVGLLAYVLVLGPGPPPLGPLARALQLALAAFQLLNLLGNVG LFLRSDPSIRGVMLAGRGLGQGWAYCYQCQSQVPPRSGHCSACRVCILRRDHHCRLLGRCVGFGNYPFL CLLLHAAGVLLHVSVLLGPALSALLRAHTPLHMAALLLLPWLMLLTGRVSLAQFALAFVTDTCVAGALLC GAGLLFHGMLLLRGQTTWEWARGQHSYDLGPCHNLQAALGPRWALVWLWPFLASPLPGDGITFQTTADVG HTAS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_207340
ORF Size:	852 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_207340.1
RefSeq Size:	1319 bp
RefSeq ORF:	855 bp
Locus ID:	254359
UniProt ID:	Q6UX98
Cytogenetics:	11q13.2
MW:	30.2 kDa