

Product datasheet for **RC209320L3V**

TRA16 (NR2C2AP) (NM_176880) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TRA16 (NR2C2AP) (NM_176880) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRA16
Synonyms:	TRA16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_176880
ORF Size:	417 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209320).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_176880.4
RefSeq Size:	1322 bp
RefSeq ORF:	420 bp
Locus ID:	126382
UniProt ID:	Q86WQ0
Cytogenetics:	19p13.11
Protein Families:	Druggable Genome
MW:	15.9 kDa


[View online »](#)

Gene Summary:

May act as a repressor of NR2C2-mediated transactivation by suppressing the binding between NR2C2/TR4 and the TR4-response element in target genes.[UniProtKB/Swiss-Prot Function]