

Product datasheet for **RC209323L3V**

FAM171B (NM_177454) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM171B (NM_177454) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM171B
Synonyms:	KIAA1946
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177454
ORF Size:	2484 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209323).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177454.3 , NP_803237.2
RefSeq Size:	3557 bp
RefSeq ORF:	2481 bp
Locus ID:	165215
UniProt ID:	Q6P995
Cytogenetics:	2q32.1
Protein Families:	Transmembrane
MW:	92.5 kDa



[View online »](#)