

## Product datasheet for **RC209344A1V**

### Human CCL28 (NM\_148672) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human CCL28 (NM\_148672) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** CCL28  
**Synonyms:** CCK1; MEC; SCYA28  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC209344 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCAGCAGAGAGGACTCGCCATCGTGGCCTTGCTGTCTGTGCGGCCCTACATGCCTCAGAAGCCATAC  
TTCCATTGCCTCCAGCTGTTGCACGGAGGTTTCACATCATATTTCCAGAAGGCTCCTGGAAAGAGTGAA  
TATGTGTGCGATCCAGAGAGCTGATGGGGATTGTGACTTGGCTGCTGTATCCTTCATGTCAAGCGCAGA  
AGAATCTGTGTGAGCCCGCACACCATACTGTTAAGCAGTGGATGAAAGTCAAGCTGCCAAGAAAAATG  
GTAAAGGAAATGTTTGCCACAGGAAGAAACACCATGGCAAGAGGAACAGTAACAGGGCACATCAGGGGAA  
ACACGAAACATACGGCCATAAACTCCTTAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC209344 protein sequence  
Red=Cloning site Green=Tags(s)

MQQRGLAIVALAVCAALHASEAILPIASSCCTEVSHHISRRLRERVMCRIQRADGDCDLAAVILHVKRR  
RICVSPHNHTVKQWMKVQAANKNGKGNVCHRKKHHGKRNSNRAHQGKHETYGHKTPY

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_148672



[View online »](#)

<b>ORF Size:</b>	381 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_148672.2</a></u> , <u><a href="#">NP_683513.1</a></u>
<b>RefSeq Size:</b>	3126 bp
<b>RefSeq ORF:</b>	384 bp
<b>Locus ID:</b>	56477
<b>UniProt ID:</b>	<u><a href="#">Q9NRJ3</a></u>
<b>Cytogenetics:</b>	5p12
<b>MW:</b>	14.3 kDa