

Product datasheet for **RC209347L3V**

PFDN5 (NM_002624) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PFDN5 (NM_002624) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PFDN5
Synonyms:	MM-1; MM1; PFD5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002624
ORF Size:	462 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209347).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002624.2
RefSeq Size:	743 bp
RefSeq ORF:	465 bp
Locus ID:	5204
UniProt ID:	Q99471
Cytogenetics:	12q13.13
Domains:	DUF232
Protein Families:	Transcription Factors



[View online »](#)

MW: 17.3 kDa

Gene Summary: This gene encodes a member of the prefoldin alpha subunit family. The encoded protein is one of six subunits of prefoldin, a molecular chaperone complex that binds and stabilizes newly synthesized polypeptides, thereby allowing them to fold correctly. The complex, consisting of two alpha and four beta subunits, forms a double beta barrel assembly with six protruding coiled-coils. The encoded protein may also repress the transcriptional activity of the proto-oncogene c-Myc. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]