

Product datasheet for **RC209368L4V**

ARL8B (NM_018184) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ARL8B (NM_018184) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARL8B
Synonyms:	ARL10C; Gie1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018184
ORF Size:	558 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209368).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018184.2
RefSeq Size:	3008 bp
RefSeq ORF:	561 bp
Locus ID:	55207
UniProt ID:	Q9NVJ2
Cytogenetics:	3p26.1
Domains:	RAB, SAR, ARF, arf
MW:	21.5 kDa



[View online »](#)

Gene Summary:

Plays a role in lysosome motility (PubMed:16537643, PubMed:25898167). In neurons, mediates the anterograde axonal long-range transport of presynaptic lysosome-related vesicles required for presynaptic biogenesis and synaptic function (By similarity). May play a role in chromosome segregation (PubMed:15331635).[UniProtKB/Swiss-Prot Function]