

Product datasheet for **RC209369L3V**

TIMM23 (NM_006327) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TIMM23 (NM_006327) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TIMM23
Synonyms:	TIM23
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006327
ORF Size:	1449 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209369).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006327.2
RefSeq Size:	1239 bp
RefSeq ORF:	630 bp
Locus ID:	100287932
UniProt ID:	O14925
Cytogenetics:	10q11.22
MW:	21.8 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is part of a complex located in the inner mitochondrial membrane that mediates the transport of transit peptide-containing proteins across the membrane. Multiple transcript variants, one protein-coding and others not protein-coding, have been found for this gene. [provided by RefSeq, Jul 2012]