

Product datasheet for **RC209373A1V**

Human RPS17 (NM_001021) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RPS17 (NM_001021) AAV Particle
Tag: Myc-DDK
Symbol: RPS17
Synonyms: DBA4; RPS17L; RPS17L1; RPS17L2; S17
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC209373 representing NM_001021
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGCCCGTTCGCACCAAAACCGTGAAGAAGCGGCCCGGGTCATCATAGAAAAGTACTACACGCGCC
TGGGCAACGACTTCCACACGAACAAGCGCGTGTGCGAGGAGATCGCCATTATCCCAGCAAAAAGCTCCG
CAACAAGATAGCAGGTTATGTCACGCATCTGATGAAGCGAATTCAGAGAGGCCAGTAAGAGGTATCTCC
ATCAAGCTGCAGGAGGAGAGAGAGAAAAGGAGAGACAATTATGTTCTGAGGTCTCAGCCTTGGATCAGG
AGATTATTGAAGTAGATCCTGACACTAAGGAAATGCTGAAGCTTTGGACTTCGGCAGTCTGTCCAACCT
TCAGGTCACCTCAGCCTACAGTTGGGATGAATTTCAAAACGCCTCGGGGACCTGTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC209373 representing NM_001021
Red=Cloning site Green=Tags(s)

MGRVRTKTVKKAARVIEKYYTRLGNDFHTNKRVCEEIAIIPSKLRNKIAGYVTHLMKRIQRGPVIRGIS
IKLQEEERERRDNYVPEVSALDQEIIEVDPDTKEMLKLLDFGSLSNLQVTQPTVGMNFKTPRGPV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001021



[View online »](#)

ORF Size:	405 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001021.3</u>
RefSeq Size:	590 bp
RefSeq ORF:	408 bp
Locus ID:	6218
UniProt ID:	<u>P08708</u>
Cytogenetics:	15q25.2
MW:	15.6 kDa