

## Product datasheet for **RC209391L1V**

### **KK LC 1 (CT83) (NM\_001017978) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KK LC 1 (CT83) (NM_001017978) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KK LC 1
Synonyms:	CXorf61; KK-LC-1; KKLC1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001017978
ORF Size:	339 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209391).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001017978.1</a>
RefSeq Size:	580 bp
RefSeq ORF:	342 bp
Locus ID:	203413
UniProt ID:	<a href="#">Q5H943</a>
Cytogenetics:	Xq23
Protein Families:	Transmembrane
MW:	12.8 kDa


[View online »](#)