

Product datasheet for **RC209395A1V**

Human RHOJ (NM_020663) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RHOJ (NM_020663) AAV Particle
Tag: Myc-DDK
Symbol: RHOJ
Synonyms: ARHJ; RASL7B; TC10B; TCL
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC209395 ORF sequence
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAAGTCAAAGAGGGAAGTACAGCAGCTGCGGCTGCAGGGGCAACGACGAGAAGAAGATGTTGAAGT
GTGTGGTGGTGGGGACGGTGCCGTGGGAAAACCTGCCTGCTGATGAGCTACGCCAACGACGCCTTCCC
AGAGGAATACGTGCCACTGTGTTTGACCACTATGCAGTTACTGTGACTGTGGGAGGCAAGCAACACTTG
CTCGGACTGTATGACACCGCGGGACAGGAGGACTACAACAGCTGAGGCCACTCTCTACCCCAACACGG
ATGTGTTTTTGATCTGCTTCTCTGCTAAACCTGCCTCTTACCACAATGTCCAGGAGGAATGGGTCCC
CGAGCTCAAGGACTGCATGCCTCACGTGCCTTATGTCCTCATAGGGACCCAGATTGATCTCCGTGATGAC
CCAAAAACCTTGCCCGTTTGCTGTATATGAAAGAGAAACCTCTCACTTACGAGCATGGTGTGAAGCTCG
CAAAAGCGATCGGAGCACAGTGTACTTGAATGTTTCAGCTCTGACTCAGAAAGGTCTCAAAGCGGTTTT
TGATGAAGCAATCCTCACCATTTCACCCCAAGAAAAAGAAACGCTGTTCTGAGGGTCACAGCTGC
TGTTCAATTATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC209395 protein sequence
Red=Cloning site **Green**=Tags(s)

MNCKEGTDSSCGRCRGNDEKKMLKCVVVG DGAVGKTCLLSYANDAFPEEYVPTVFDHYAVTVTVGGKQHL
LGLYDTAGQEDYNQLRPLSYNTDVFLICFSVNPASYHNVQE EWVPELKDCMPHVYPYVLIGTQIDLRDD
PKTLARLLYMKEKPLTYEHGVLAKAIGAQCYLECSALTQKGLKAVFDEAILTFHPKKKKKRCSGHS
CSII

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_020663
ORF Size:	642 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_020663.2
RefSeq Size:	3619 bp
RefSeq ORF:	645 bp
Locus ID:	57381
UniProt ID:	Q9H4E5
Cytogenetics:	14q23.2
MW:	23.8 kDa