

Product datasheet for **RC209398A1V**

Human CRSP8 (MED27) (NM_004269) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CRSP8 (MED27) (NM_004269) AAV Particle
Tag:	Myc-DDK
Symbol:	CRSP8
Synonyms:	CRAP34; CRSP8; CRSP34; MED3; NEDSCAC; TRAP37
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC209398 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGACGTGATAAATGTCAGTGTGAACCTGGAGGCCTTTTCCCAGGCCATTAGTGCCATCCAGGCGC
TGCGCTCCAGCGTGAGCAGGGTGTTCGACTGCCTGAAGGATGGGATGCGGAACAAGGAGACGCTGGAGGG
CCGGGAGAAGGCCTTTATTGCGCACTTCCAGGACAATTACATTCCGGTCAACCGGACCTCAATGAGCTG
GAACGCTCTGAGCAATCTGGTAGGCAAGCCATCTGAGAACCATCCTCTTCATAACAGTGGGCTGTTAAGCC
TGGATCCTGTGCAGGACAAAACCTCTCTATAGTCAACTCCTTCAAGCATATAAGTGGTCAAACAAGTT
GCAGTACCATGCAGGACTAGCATCTGGCCTTTTAAATCAGCAGTCATTGAAGCGTTCCGCTAATCAGATG
GGAGTATCTGCCAAACGTAGACCAAAGGCTCAGCCACAACTCTTGCTTACCACCTCAATATGTTGATG
ATGTGATCAGCCGATTGACAGGATGTTTCTGAAATGTCCATCCACTTATCCAGACCAATGGAACATC
AGCAATGCTTCTGGTGACCTTGGGAAAGGTGTTGAAAGTGATCGTCGTCATGCGGAGCCTGTTTATTGAT
CGAACAATAGTAAAGGGATATAACGAGAATGTCTACACAGAAGATGGCAAGCTTGATATATGGTCCAAAT
CCAATATCAAGTATTCCAGAAGGTGACAGACCATGCCACCACTGCCCTGCTCCACTATCAGTGCCCCA
GATGCCGATGTCGTGGTCCGATCCTTCATGACCTGGTTAAGAAGTTACATAAAGCTGTTCCAGGCCCG
TGCCAGCGCTGCGGGAAGTTTCTGCAGGACGGCCTTCCCCGACATGGAGGGATTCCGAACCCTCGAAG
CCTTCCATGACACCTGCCGCGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC209398 protein sequence Red=Cloning site Green=Tags(s) MADVINSVNLEAFSQAISAIQALRSSVSRVFDCLKDGMRNKETLEGREKAFIAHFQDNLHNSVNRDLNEL ERLSNLVGKPSENHPLHNSGLLSLDPVQDKTPLYSQLQAYKWSNKLQYHAGLASGLLNQQSLKRSANQM GVSARRPKAQPTTLVLPQYVDDVISRIDRMFPMSIHL SRPNGTSAMLLVTLGKVLKVIIVMRS LFID RTIVKGYNENVYTEDGKLDIWSKSNYQVFQKVTDHATTALLHYQLPQMPDVVVR SFMTWLR SYIKLFQAP CQRCGKFLQDGLPPTWRDFRTLEAFHDTCRQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_004269
ORF Size:	933 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_004269.2</u>
RefSeq Size:	1418 bp
RefSeq ORF:	936 bp
Locus ID:	9442
UniProt ID:	<u>Q6P2C8</u>
Cytogenetics:	9q34.13
MW:	35.4 kDa